




## Service manual for VELA Person Lift

Identification number: \_\_\_\_\_

## Checklist:

1. Check the lift for cracks, fractures and deformation.
2. Check that the person lift can lift a person and a body support of 125 kg in total. One complete cycle must be run with maximum load. The person lift is adjusted all the way up and down again. Listen for irregular noises during use.
3. Check that the emergency stop works (A).
4. Retighten all bolts (B), check for backlash and see that joints are running smoothly.
5. Grease the moving joints and retighten the bolts (C).
6. Check bearing bushings for wear (D).
7. Check the power cord for wear.
8. Check that the charger is working by connecting it and see that it charges.



Inspection notes: \_\_\_\_\_

\_\_\_\_\_

Other comments: \_\_\_\_\_

\_\_\_\_\_

Inspection date: \_\_\_\_\_

Next inspection date: \_\_\_\_\_

Carried out by: \_\_\_\_\_

Copy for the VELA Person Lift user